



MRI READY LEADLESS SYSTEM:

MRI Scan Clinician Checklist



PATIENT: _____ DATE OF BIRTH: _____

DEVICE INFORMATION

PRODUCT: _____

MANUFACTURER: _____ MODEL #: _____

STEPS FOR AN MRI SCAN¹ CARDIOLOGY DEPARTMENT

Please refer to the Abbott MRI Ready Leadless System Manual at <https://medical.abbott/manuals>

MRI SYSTEM CONDITIONS	YES	NO
Patient has an MR Conditional Aveir™ Leadless Pacemaker.		
Confirmed that no adverse conditions to MRI scanning are present.		
Confirmed the absence of additional cardiac hardware.		
MR Conditional Aveir Leadless Pacemaker is implanted in the right ventricle.		
Leadless pacemaker has sufficient battery (not at end-of-life EOL).		
PRE-MRI SCAN DEVICE SETTINGS	YES	NO
Confirmed pacing capture thresholds are stable.		
Confirmed pacing capture threshold ≤ 2.5 V at a pulse width of 0.5 ms.		
Confirmed MRI Settings status and the programmed MRI Mode setting. Please select a pacing rate to avoid competitive pacing: VOO (Bipolar) Pacing rate _____ bpm Pacing Off _____		
<i>REMINDER: It is required to program the device to MRI Settings as part of the MRI scan workflow.</i>		
Confirmed no complaints of diaphragmatic stimulation at a pacing output of 5.0 V at a pulse width of 1.0 ms if the patient is programmed to an asynchronous pacing mode when MRI Settings are enabled.		
Printed the MRI Summary Report of the patient's MRI Settings and permanently programmed parameters.		

STEPS FOR AN MRI SCAN[†] CARDIOLOGY DEPARTMENT

Please refer to the Abbott MRI Ready Leadless System Manual at <https://medical.abbott/manuals>

POST-MRI SCAN DEVICE SETTINGS	YES	NO
Disabled the MRI Settings on the leadless pacemaker, restoring the permanently programmed settings.		
Confirmed the permanently programmed settings are appropriate.		
Checked the pacing capture thresholds after the scan to ensure that the pacing parameters are programmed adequately for the patient based on the threshold.		

PHYSICIAN NAME: _____

DATE: _____

PHYSICIAN SIGNATURE: _____

STEPS FOR AN MRI SCAN[†] RADIOLOGY DEPARTMENT

Please refer to the Abbott MRI Ready Leadless System Manual at <https://medical.abbott/manuals>

MRI SCAN CONDITIONS	YES	NO
Confirmed patient has an MR Conditional Aveir Leadless Pacemaker.		
Confirmed that no adverse conditions to MRI scanning are present.		
Selected the correct Scan Parameters for your MRI equipment.		
Confirmed the MRI Settings status (MRI Settings must be enabled before the scan).		
Performed the scan and monitor the patient.		
After patient receives MRI scan, the leadless pacemaker must be interrogated and programmed by the patient's device management physician or clinician to disable MRI Settings.		

NAME: _____

DATE: _____

SIGNATURE: _____

TECHNICAL SUPPORT

Abbott maintains 24-hour phone lines for technical questions and support: 1-818-362-6822

1-800-722-3774 (toll-free within North America)

+46 8 474 4147 (Sweden)

+ 61 2 9936 1200 (Australia)

medical.abbott/manuals

For additional assistance, call your local Abbott representative.

REFERENCES

- a. For any additional information, including warnings, precautions, and potential adverse events, please refer to the Abbott MRI Ready Leadless System Manual at <https://medical.abbott/manuals> or visit cardiovascular.abbott/mriready.
- b. Do not bring any external devices, such as a Merlin™ PCS Programmer, Aveir Link Module, or a smart phone or tablet running a mobile application, into the scanner magnet room (Zone IV).

This document is intended to be used as a guide for health care facilities.

Abbott

15900 Valley View Court Sylmar, CA 91342

Tel: +1 818 362 6822

[Abbott.com](https://www.abbott.com)

Rx Only

Brief Summary: Prior to using these devices, please review the Instructions for Use for a complete listing of indications, contraindications, warnings, precautions, potential adverse events and directions for use.

™ Indicates a trademark of the Abbott group of companies.

© 2022 Abbott. All Rights Reserved.
MAT-2005394 v2.0 | Item approved for U.S. use.

