




**CHANGE THE COURSE OF HEART FAILURE
FOR SELECT HEART FAILURE PATIENTS WITH CLINICALLY
SIGNIFICANT SECONDARY MITRAL REGURGITATION**

**TRANSCATHETER MITRAL VALVE REPAIR
WITH MITRACLIP™ THERAPY**



VASTLY UNDERTREATED. IMPACTING SURVIVAL.


HIGHLY PREVALENT¹



1 in 5 of your heart failure patients have moderate-to-severe and severe secondary MR.^{1-6*}

VASTLY UNDERTREATED

SURGICAL INTERVENTION NOT OFFERED OR DENIED**




49% of patients with symptomatic severe MR were not operated due to age, co-morbidities, or impaired LV.⁷

SIGNIFICANT IMPACT ON MORTALITY



If left untreated, MR initiates a cascade of events leading to death.

SEVERE SECONDARY MR IS AN INDEPENDENT PREDICTOR OF MORTALITY¹¹



Heart failure (HF) patients with severe secondary MR have an even poorer prognosis.¹¹

WHAT ARE THE SYMPTOMS OF MR?

In some cases, patients may have MR but not experience any symptoms. In other cases, patients may develop symptoms of heart failure such as:¹²⁻¹⁴



- Fatigue
- Inability to exercise
- Decrease in appetite
- Dry, hacking cough (often worse when lying down)
- Shortness of breath (especially at night)
- Fainting
- Weight gain from retaining fluid
- Accumulation of fluid in feet, ankles, and lungs (edema)

WHAT ARE THE TREATMENT OPTIONS?



MEDICATIONS

A number of medical therapies are available that can treat the symptoms of MR, but they do not address the underlying pathology.¹⁵ Patients with heart failure and clinically significant mitral regurgitation are on guideline-directed medical therapy (GDMT).



SURGERY (PRIMARY MR)

Surgical management, in the form of mitral valve repair or replacement, is the standard of care for patients with severe primary mitral regurgitation. However, up to 50% of patients with MR may not meet the eligibility criteria due to risks associated with age or the presence of comorbidities.¹⁶ For patients with severe secondary mitral regurgitation, guideline-directed medical therapy is the current standard of care.



CARDIAC RESYNCHRONISATION THERAPY (SECONDARY MR)

Cardiac resynchronisation therapy is an option recommended for select HF patients.⁷




TRANSCATHETER MITRAL VALVE REPAIR (TMVr) WITH MITRACLIP

TMVr is a minimally invasive mitral valve repair option for patients with MR who are not considered to be suitable candidates for surgery.¹² MitraClip is the only minimally invasive repair option that is included in the ESC heart valve guidelines.¹²⁻¹⁴

*Heart Failure patients with reduced EF and with moderate to severe and severe secondary MR.
**Based on a survey of severely symptomatic MR patients in NYHA Class III-IV (n = 396); 10% had surgery the following year. The remainder had no surgery; medical management only.

A LANDMARK STUDY. TO ESTABLISH A NEW STANDARD OF CARE.¹⁷



Published in the New England Journal of Medicine, the COAPT Trial was a parallel-controlled, open-label, multicenter trial in 614 patients with heart failure and moderate-to-severe (3+) or severe (4+) secondary MR who remained symptomatic despite maximally-tolerated GDMT.¹⁷


PRIMARY ENDPOINTS:

Effectiveness: All HF hospitalizations through 24 months, analyzed when last patient completes 12-month follow-up

Reduces Complications: 96.6% reported freedom from device-related complications through 12 months


MR KILLS—THE COAPT TRIAL SHOWS THAT IT IS NOT JUST A MARKER

MITRACLIP THERAPY SAVES LIVES¹⁸



33%
RELATIVE RISK REDUCTION IN MORTALITY

NUMBER NEEDED TO TREAT (NNT) TO PREVENT ONE DEATH FROM ANY CAUSE*




7.9
NUMBER NEEDED TO TREAT TO PREVENT ONE DEATH

*Data from different trials with similar follow up periods; incremental benefits due to test drug/device above background therapy.
**Median follow-up duration.

MITRACLIP IS THE ONLY DEVICE THERAPY SHOWN TO IMPROVE SURVIVAL OF HEART FAILURE PATIENTS WITH SECONDARY MR¹⁷


REDUCES HOSPITALIZATIONS FOR HF¹⁸



51%
RELATIVE RISK REDUCTION IN HEART FAILURE HOSPITALIZATIONS

3.0
NUMBER NEEDED TO TREAT TO PREVENT ONE HEART FAILURE HOSPITALIZATION


IMPROVES QUALITY OF LIFE²¹



2.5X
MORE LIKELY TO EXPERIENCE A LARGE IMPROVEMENT IN QUALITY OF LIFE WITH MITRACLIP

Note: KCCQ Minimum for Clinically Important Difference (MCID)= 5 points; Large Improvement Defined as ≥ 20 Points in KCCQ from baseline; quality of life is assessed only in surviving patients

PROVIDES DURABLE MR REDUCTION¹⁷



99.1%
OF MITRACLIP PATIENTS HAD MR $\leq 2+$ AT 24 MONTHS


Note: Unpaired data

96.6%

FREEDOM FROM DEVICE RELATED COMPLICATIONS AT 12 MONTHS

WHAT IS THE MITRA-FR TRIAL AND WHY ARE THE RESULTS DIFFERENT FROM COAPT?

THE MITRA-FR TRIAL²²



- The MITRA-FR trial was an independent multicenter, randomized, open-label, controlled phase 3 trial that was conducted in France
- Hospices Civils de Lyon, a public academic institution assumed overall responsibility for the trial
- Primary funding was provided by the French Ministry of Health and Research National Program
- Abbott provided the devices

PRIMARY ENDPOINT:

- All cause deaths or unplanned hospitalization for heart failure at 12 month


MITRA-FR AND COAPT ARE 2 DIFFERENT TRIALS THAT STUDIED 2 DIFFERENT PATIENT GROUPS²³

| | MITRA-FR | COAPT |
|----------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------|
| Severe MR entry criteria | Severe SMR by EU guidelines: EROA >20 mm ² or RV >30 mL/beat | Severe SMR by US guidelines: EROA >30 mm ² or RV >45 mL/beat |
| % of patients with an EROA <0.3 cm² | 52% | 14% |
| LVEDV | 70% of patients had an LV EDD >65 mm | Exclusion criteria of patients with LVESD >70mm |
| GDMT at baseline and FU | Receiving HF meds at baseline—allowed variable adjustment in each group during follow-up per “real-world” practice | CEC confirmed pts were failing maximally-tolerated GDMT at baseline—few major changes during follow-up |
| Acute results: No clip / ≥3+ MR | 9% / 9% | 5% / 5% |
| Procedural complications* | 14.6% | 8.5% |
| 12-mo MitraClip ≥3+ MR | 17% | 5% |
| Mortality outcomes | No difference in device vs control group | 38% reduction |

*MITRA-FR definition: device implant failure, transfusion or vascular complications requiring surgery, ASD, card shock, cardiac embolism/stroke, tamponade, urgent cardiac surgery

HOWEVER, QUALITY OF LIFE IMPROVED IN THE MITRACLIP GROUP²²

The Quality of Life score (according to the EQ5D scale developed by the EuroQol group) and the 6 Minute Walk Test Distance were higher at 12 months in the MitraClip group than in the control group at 12 months.




CLINICAL BENEFITS FOR BROAD RANGE OF PATIENTS WITH SMR ARE DEMONSTRATED WITH MITRACLIP²²

| PATIENT SUBGROUP | | Clinical Outcomes |
|-----------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------|
| EROA > 0.30 cm² or LVEDVi ≤ 96 mL/m² [90% of COAPT Pts] | COAPT ¹⁷ | Reduces MR, reduces heart failure hospitalizations, improves survival, quality of life, and functional capacity |
| EROA ≤ 0.30 cm² or LVEDVi > 96 mL/m² [10% of COAPT Pts] | COAPT ¹⁷ (MITRA-FR like) | Improves quality of life and functional capacity |
| EROA > 0.20 cm² LVEDVi 135 ± 35 mL/m² | MITRA-FR ²² | Reduces MR, improves quality of life and functional capacity* |
| Broad range of patients with SMR | EVEREST II ²⁴⁻²⁵ , ACCESS-EU ²⁶⁻²⁸ , REALISM ²⁹⁻³⁰ , TRAMI ³¹⁻³³ | Reduces MR, reduces heart failure hospitalizations, improves quality of life, and functional capacity |

*Table 56, Supplement to Obadia JF et al. Percutaneous repair or medical treatment for secondary regurgitation. *N Engl J Med* 2018; 379:2297-2306.

GLOBALLY, MITRACLIP HAS OVER 16 YEARS OF CLINICAL EXPERIENCE³⁴



TAKING ACTION FOR MR PATIENTS

PATIENT TREATMENT PATHWAY

1 IDENTIFY MITRACLIP CANDIDATES FOR EVALUATION

Candidates for MitraClip therapy should meet the following criteria:³⁵

- ✓ Significant MR
- ✓ Severely symptomatic
- ✓ Symptoms and MR severity persist despite maximally-tolerated GDMT
- ✓ Not suitable for surgery because of age or severe comorbidities such as liver disease and pulmonary hypertension



A diagnosis of severe MR is based on multiple parameters, including:¹⁹

- Valve anatomy
- Regurgitant severity
- Left ventricular size
- Systolic function


2 SEND TRANSTHORACIC ECHO (TTE) TO ASSIST IN EVALUATIONS

Send patients to a MitraClip Center with TTE views to ensure efficient analysis of the mitral valve and to assess initial anatomic eligibility of MitraClip therapy.



3 CONSULT A TRANSCATHETER MITRAL VALVE REPAIR (TMV_r) CENTER

The 2016 ESC/HF and 2017 ESC/EACTS guidelines recommend that complex patients like those with severe, symptomatic MR be evaluated by a Heart Valve Center of Excellence, such as a MitraClip Center.¹⁹



ASK YOUR MITRACLIP REPRESENTATIVE FOR A LIST OF TMV_r CENTERS NEAR YOU.

ABBOTT PARTNERS WITH YOU TO PROVIDE THE MOST COMPREHENSIVE SET OF SOLUTIONS FOR PATIENTS WITH MR

WITH SUPPORT FROM ABBOTT STRUCTURAL HEART, YOU CAN:



Identify Patients Efficiently



Provide Patient Education



Connect with TMVr Centers

If your patient is a candidate for MitraClip, refer them to a specialist for treatment, or visit AdvancedHeartTherapies.com for more information.

References

1. AHA Heart Disease and Stroke Statistics Update, *Circulation* 2017. 2. Yancy CW et al, *JACC* 2013. 3. Pecini et al, *EHJ* 2011. 4. Asgar et al, *JACC* 2015. 5. Nieminen et al, *EHJ* 2006. 6. Patel et al, *Journal of Cardiac Failure* 2004. 7. Mirabel M, Iung B, Baron G, et al. What are the characteristics of patients with severe, symptomatic, mitral regurgitation who are denied surgery? *Eur Heart J*. 2007;28(11):1358-1365. 8. Enriquez-Sarano M, Avierinos JF, Messika-Zeitoun D, et al. Quantitative determinants of the outcome of asymptomatic mitral regurgitation. *N Engl J Med*. 2005;352(9):875-883. 9. Cioffi G, Tarantini L, De Feo S, et al. Functional mitral regurgitation predicts 1-year mortality in elderly patients with systolic chronic heart failure. *Eur J Heart Fail*. 2005;7(7):1112-1117. 10. Grigioni F, Tribouilloy C, Avierinos JF, et al; MIDA Investigators. Outcomes in mitral regurgitation due to flail leaflets: a multicenter European study. *JACC Cardiovasc Imaging*. 2008;1(2):133-141. 11. Goliash G, et al. Refining the prognostic impact of functional mitral regurgitation in chronic heart failure. *Eur Heart J*. 2018 Jan 1;39(1):39-46. 12. Nishimura RA, Otto CM, Bonow RO, et al. 2017 AHA/ACC focused update of the 2014 AHA/ACC guideline for the management of patients with valvular heart disease: a report of the American College of Cardiology/American Heart Association Task Force on Clinical Practice Guidelines. *Circulation*. 2017;136(9):1-123. DOI: 10.1161/CIR.0000000000000503. 13. Baumgartner H, Falk V, Bax JJ, et al. 2017 ESC/EACTS Guidelines for the management of valvular heart disease: The Task Force for the Management of Valvular Heart Disease of the European Society of Cardiology (ESC) and the European Association for Cardio-Thoracic Surgery (EACTS). *Eur Heart J*. 2017;00:1-53. 14. Ponikowski P, Voors AA, Anker SD, et al. 2016 ESC Guidelines for the diagnosis and treatment of acute and chronic heart failure: The Task Force for the Diagnosis and Treatment of Acute and Chronic Heart Failure of the European Society of Cardiology (ESC). Developed with the special contribution of the Heart Failure Association (HFA) of the ESC. *Eur J Heart Fail*. 2016;18(8):891-975. 15. A. Young and T. Feldman. *Current Cardiology Reports*. 16, 443 (2014). 16. Joint Task Force on the Management of Valvular Heart Disease of the European Society of Cardiology, et al. *European Heart Journal*. 2012;33:2451. 17. Stone GW, Lindenfeld J, Abraham WT, et al. Transcatheter mitral-valve repair in patients with heart failure. *N Engl J Med*. September 23, 2018. 18. Mack M. COAPT: Three-year outcomes from a randomized trial of transcatheter mitral valve leaflet approximation in patients with heart failure and secondary mitral regurgitation. Presented at TCT 2019. 19. McMurray JVV, Packer M, Desai AS, et al. *N Engl J Med* 2014;371:993-1004. 20. Manually calculated from Table 2: Cleland JG et al. *N Engl J Med* 2005;352:1539-1549. 21. Arnold SV et al. Health status after transcatheter mitral valve repair in heart failure and secondary mitral regurgitation. *JACC* Mar 2019, 25951; DOI: 10.1016/j.jacc.2019.02.010. 22. Obadia JF et al. Percutaneous repair or medical treatment for secondary mitral regurgitation. *N Engl J Med* 2018; 379:2297-2306. 23. Stone G. Pivotal transcatheter FMR device trials: Focus on COAPT and MITRA-FR, with implications for other transcatheter mitral valve device investigations. FDA Townhall, TCT 2018. 24. Glower DD, Kar S, Trento A, Lim DS, Bajwa T, Quesada R, et al. Percutaneous mitral valve repair for mitral regurgitation in high risk patients. *JACC*. 2014;64(2):172-81. 25. Whitlow PL, Feldman T, Pedersen WR, Lim DS, Kipperman R, Smalling R, et al. Acute and 12-month results with catheter based mitral valve leaflet repair: the EVEREST II (Endovascular Valve Edge-to-Edge Repair) High Risk Study. *JACC*. 2012 Jan 10;59(2):130-9. PubMed PMID: 22222076. 26. Maisano F, Franzen O, Baldus S, et al. Percutaneous mitral valve interventions in the real world: early and 1-year results from the ACCESS-EU, a prospective, multicentre, nonrandomized post-approval study of the MitraClip Therapy in Europe. *JACC*. 2013;62(12):1052-61. 27. Butter C, Franzen O, Treede H, et al. Analysis of outcomes at 12 months by mitral regurgitation etiology in the ACCESS-EUROPE phase I study. In: EuroIntervention, editor. EuroPCR; Paris: EuroIntervention; 2013. 28. Hoffman I. Analysis of outcomes at 12 months by gender in the ACCESS-EUROPE Phase I Study. Presented at PCR London Valves 2013, London, UK. 2013. 29. Kar S, Lim D, Smalling R, Whisenant B, Rammohan C, Fail PS, et al. The EVEREST II REALISM continued access study: effectiveness of transcatheter reduction of significant mitral regurgitation in surgical candidates. *JACC*. 2013;61(10):E1959. 30. Lim DS, Reynolds MR, Feldman T, Kar S, Herrmann HC, Wang A, et al. Improved functional status and quality of life in prohibitive surgical risk patients with degenerative mitral regurgitation following transcatheter mitral valve repair with the MitraClip System. *JACC*. 2013 Oct 24. PubMed PMID: 24184254. 31. Puls M, Lubos E, Boekstegers P, et al. One-year outcomes and predictors of mortality after MitraClip therapy in contemporary clinical practice: results from the German transcatheter mitral valve interventions registry. *European Heart Journal*. 2016;37(8):703-12. 32. Wiebe J, Franke J, Lubos E, Boekstegers P, Schillinger W, Ouarrak T, et al. Percutaneous mitral valve repair with the MitraClip system according to the predicted risk by the logistic EuroSCORE: Preliminary results from the German Transcatheter Mitral Valve Interventions (TRAMI) registry. *Catheterization and Cardiovascular Interventions*. 2014 Mar 24. PubMed PMID: 24664460. 33. Rudolph V, Huntgeburth M, Von Bardeleben RS, et al. Clinical outcome of critically ill, not fully recompensated, patients undergoing MitraClip therapy. *European Journal Heart Failure*. 2014;16:1223-9. 34. Abbott data on file as of Nov 2019. 35. MitraClip NT XTR Clip Delivery System US Instructions for Use.

CAUTION: This product is intended for use by or under the direction of a physician. Prior to use, reference the Instructions for Use provided inside the product carton (when available) or at eifu.abbottvascular.com or at medical.abbott/manuals for more detailed information on Indications, Contraindications, Warnings, Precautions and Adverse Events. Information contained herein for DISTRIBUTION in Australia and New Zealand ONLY. Check the regulatory status of the device in areas where CE marking is not the regulation in force.

Illustrations are artist's representations only and should not be considered as engineering drawings or photographs. Photos on file at Abbott.

Abbott Australasia Pty Ltd: 299 Lane Cove Road, Macquarie Park, NSW, Australia 2113, Tel: 1800 550 939
Abbott Laboratories (NZ) Limited: Building D, 4 Pacific Drive, Mount Wellington, Auckland 1641, Tel: 0800 827 285

™ indicates a trademark of the Abbott Group of Companies.
www.Cardiovascular.Abbott

©2020 Abbott. All rights reserved. ANZ SH0144-EN 01/20 | Item approved for use in ANZ only.

