



# A LANDMARK STUDY TO ESTABLISH A NEW STANDARD OF CARE


TRANSCATHETER MITRAL VALVE REPAIR  
WITH MITRACLIP™ THERAPY



MITRAL REGURGITATION:

# VASTLY UNDERTREATED. IMPACTING SURVIVAL.

## HIGHLY PREVALENT<sup>1</sup>



The prevalence of mitral valve disease is **2-3x** greater than that of aortic valve disease, with over 4 million patients suffering from moderate or severe MR in the U.S. alone.<sup>1-3</sup>

## SIGNIFICANT IMPACT ON MORTALITY




If left untreated, MR initiates a cascade of events leading to death.

<sup>1</sup>Based on a survey of severely symptomatic MR patients in NYHA Class III-IV (n = 396); 10% had surgery the following year. The remainder had no surgery; medical management only.

<sup>2</sup> INFORMATION CONTAINED HEREIN FOR DISTRIBUTION IN AUSTRALIA AND NEW ZEALAND ONLY. CHECK THE REGULATORY STATUS OF THE DEVICE IN AREAS WHERE CE MARKING IS NOT THE REGULATION IN FORCE.


## VASTLY UNDERTREATED

SURGICAL INTERVENTION NOT OFFERED OR DENIED\*



**49%** of patients with symptomatic severe MR were not operated due to age, co-morbidities, or impaired LV.<sup>4</sup>


SEVERE SECONDARY MR IS AN INDEPENDENT PREDICTOR OF MORTALITY<sup>8</sup>



Heart failure (HF) patients with severe secondary MR have an even poorer prognosis.<sup>8</sup>



# A LANDMARK STUDY IN THE MANAGEMENT OF HEART FAILURE



Published in the New England Journal of Medicine, the COAPT Trial was a parallel-controlled, open-label, multicenter trial in 614 patients with heart failure and moderate-to-severe (3+) or severe (4+) secondary MR who remained symptomatic despite maximally-tolerated GDMT.<sup>9</sup>


### PRIMARY ENDPOINTS

**Effectiveness:** All HF hospitalizations through 24 months, analyzed when last patient completes 12-month follow-up

**Reduces Complications:** 96.6% reported freedom from device-related complications through 12 months

## MR KILLS—THE COAPT TRIAL SHOWS THAT IT IS NOT JUST A MARKER

MITRACLIP THERAPY SAVES LIVES<sup>10</sup>



**33%**

RELATIVE RISK REDUCTION IN MORTALITY


**7.9**

NUMBER NEEDED TO TREAT TO PREVENT ONE DEATH

No. at Risk	0	6	12	18	24	30	36
MitraClip + GDMT	302	269	238	219	189	128	93
GDMT alone	312	272	223	186	145	91	70

# MITRACLIP IS THE ONLY MITRAL VALVE DEVICE SHOWN TO IMPROVE SURVIVAL OF HEART FAILURE PATIENTS WITH SECONDARY MR<sup>9</sup>


## NUMBER NEEDED TO TREAT (NNT) TO PREVENT ONE DEATH FROM ANY CAUSE\*



\*Data from different trials with similar follow up periods; incremental benefits due to test drug/device above background therapy.  
\*\*Median follow-up duration

## REDUCES HOSPITALIZATIONS FOR HF<sup>10</sup>

HEART FAILURE HOSPITALIZATIONS ARE THE LEADING COST BURDEN TO THE HOSPITAL AND HEALTH CARE SYSTEM.




**51%**  
RELATIVE RISK REDUCTION  
IN HEART FAILURE  
HOSPITALIZATIONS

**3.0**  
NUMBER NEEDED TO TREAT  
TO PREVENT ONE HEART  
FAILURE HOSPITALIZATION

\*Data from different trials with similar follow up periods; incremental benefits due to test drug/device above background therapy.


INFORMATION CONTAINED HEREIN FOR DISTRIBUTION IN AUSTRALIA AND NEW ZEALAND ONLY. CHECK THE REGULATORY STATUS OF THE DEVICE IN AREAS WHERE CE MARKING IS NOT THE REGULATION IN FORCE.

## PROVIDES DURABLE MR REDUCTION<sup>9</sup>



Note: Unpaired data

## IMPROVES QUALITY OF LIFE<sup>13</sup>




Note: KCCQ Minimum for Clinically Important Difference (MCID)= 5 points; Large Improvement Defined as ≥20 Points in KCCQ from baseline; quality of life is assessed only in surviving patients

**96.6%** FREEDOM FROM DEVICE RELATED COMPLICATIONS AT 12 MONTHS

# WHAT IS THE MITRA-FR TRIAL AND WHY ARE THE RESULTS DIFFERENT FROM COAPT?

## THE MITRA-FR TRIAL<sup>14</sup>



- The MITRA-FR trial was an independent multicenter, randomized, open-label, controlled phase 3 trial that was conducted in France
- Hospices Civils de Lyon, a public academic institution assumed overall responsibility for the trial
- Primary funding was provided by the French Ministry of Health and Research National Program
- Abbott provided the devices

### PRIMARY ENDPOINT:

- All cause deaths or unplanned hospitalization for heart failure at 12 month


## MITRA-FR AND COAPT ARE 2 DIFFERENT TRIALS THAT STUDIED 2 DIFFERENT PATIENT GROUPS<sup>15</sup>

	MITRA-FR	COAPT
<b>Severe MR entry criteria</b>	Severe SMR by EU guidelines: EROA >20 mm <sup>2</sup> or RV >30 mL/beat	Severe SMR by US guidelines: EROA >30 mm <sup>2</sup> or RV >45 mL/beat
<b>% of patients with an EROA &lt;0.3 cm<sup>2</sup></b>	52%	14%
<b>LVEDV</b>	70% of patients had an LV EDD >65 mm	Exclusion criteria of patients with LVESD >70mm
<b>GDMT at baseline and FU</b>	Receiving HF meds at baseline—allowed variable adjustment in each group during follow-up per “real-world” practice	CEC confirmed pts were failing maximally-tolerated GDMT at baseline—few major changes during follow-up
<b>Acute results: No clip / ≥3+ MR</b>	9% / 9%	5% / 5%
<b>Procedural complications*</b>	14.6%	8.5%
<b>12-mo MitraClip ≥3+ MR</b>	17%	5%
<b>Mortality outcomes</b>	No difference in device vs control group	38% reduction

\*MITRA-FR definition: device implant failure, transfusion or vascular complications requiring surgery, ASD, card shock, cardiac embolism/stroke, tamponade, urgent cardiac surgery

## HOWEVER, QUALITY OF LIFE IMPROVED IN THE MITRACLIP GROUP<sup>14</sup>

The Quality of Life score (according to the EQ5D scale developed by the EuroQol group) and the 6 Minute Walk Test Distance were higher at 12 months in the MitraClip group than in the control group at 12 months.



## CLINICAL BENEFITS FOR BROAD RANGE OF PATIENTS WITH SMR ARE DEMONSTRATED WITH MITRACLIP<sup>14</sup>

PATIENT SUBGROUP		Clinical Outcomes
<b>EROA &gt; 0.30 cm<sup>2</sup> or LVEDVi ≤ 96 mL/m<sup>2</sup> [90% of COAPT Pts]</b>	COAPT <sup>9</sup>	Reduces MR, reduces heart failure hospitalizations, improves survival, quality of life, and functional capacity
<b>EROA ≤ 0.30 cm<sup>2</sup> or LVEDVi &gt; 96 mL/m<sup>2</sup> [10% of COAPT Pts]</b>	COAPT <sup>9</sup> (MITRA-FR like)	Improves quality of life and functional capacity
<b>EROA &gt; 0.20 cm<sup>2</sup> LVEDVi 135 ± 35 mL/m<sup>2</sup></b>	MITRA-FR <sup>14</sup>	Reduces MR, improves quality of life and functional capacity*
<b>Broad range of patients with SMR</b>	EVEREST II <sup>16-17</sup> , ACCESS-EU <sup>18-20</sup> , REALISM <sup>21-22</sup> , TRAMI <sup>23-25</sup>	Reduces MR, reduces heart failure hospitalizations, improves quality of life, and functional capacity

\*Table 56, Supplement to Obadia JF et al. Percutaneous repair or medical treatment for secondary regurgitation. *N Engl J Med* 2018; 379:2297-2306.

# MITRACLIP-FIRST AND FOREMOST

## THE FIRST

- The first transcatheter mitral valve device
- Durable outcomes demonstrated to 5 years
- The first mitral valve intervention proven more effective than GDMT alone for select heart failure patients with clinically significant MR

## THE FOREMOST

- Over 16 years of dedication to the treatment of MR
- Over 100,000 patients treated worldwide
- Over 30,000 patients studied in clinical trials
- Over 16 years of dedicated support of training and education for transcatheter mitral valve procedures
- Fourth generation technology platform\*

16+

YEARS CLINICAL  
EXPERIENCE

OVER  
100K

PATIENTS TREATED  
WORLDWIDE

30K+

PATIENTS STUDIED IN  
CLINICAL TRIALS

Data on File at Abbott.

\*MitraClip G4 pending CE Mark.

For more information about Transcatheter Mitral Valve Repair with  
MitraClip Therapy or the COAPT Trial, contact your Abbott sales representative  
or visit [AdvancedHeartTherapies.com](http://AdvancedHeartTherapies.com).

## REFERENCES

1. Nkomo VT, Gardin JM, Skelton TN, Gottdiener JS, Scott CG, Enriquez-Sarano M. Burden of valvular heart diseases: a population-based study. *Lancet*. 2006;368(9540):1005-1011. 2. United States Census Bureau, 2010. Age and sex composition. Issued May 2011:1-16. 3. Abbott data on file as of Nov 2019. 4. Mirabel M, Jung B, Baron G, et al. What are the characteristics of patients with severe, symptomatic, mitral regurgitation who are denied surgery? *Eur Heart J*. 2007;28(11):1358-1365. 5. Enriquez-Sarano M, Avierinos JF, Messika-Zeitoun D, et al. Quantitative determinants of the outcome of asymptomatic mitral regurgitation. *N Engl J Med*. 2005;352(9):875-883. 6. Cioffi G, Tarantini L, De Feo S, et al. Functional mitral regurgitation predicts 1-year mortality in elderly patients with systolic chronic heart failure. *Eur J Heart Fail*. 2005;7(7):1112-1117. 7. Grigioni F, Tribouilloy C, Avierinos JF, et al; MIDA Investigators. Outcomes in mitral regurgitation due to flail leaflets: a multicenter European study. *JACC Cardiovasc Imaging*. 2008;1(2):133-141. 8. Goliash G, et al. Refining the prognostic impact of functional mitral regurgitation in chronic heart failure. *Eur Heart J*. 2018 Jan 1;39(1):39-46. 9. Stone GW, Lindenfeld J, Abraham WT, et al. Transcatheter mitral-valve repair in patients with heart failure. *N Engl J Med*. September 23, 2018. 10. Mack M. COAPT: Three-year outcomes from a randomized trial of transcatheter mitral valve leaflet approximation in patients with heart failure and secondary mitral regurgitation. Presented at TCT 2019. 11. McMurray JJV, Packer M, Desai AS, et al. *N Engl J Med* 2014;371:993-1004. 12. Manually calculated from Table 2: Cleland JG et al. *N Engl J Med* 2005;352:1539-1549. 13. Arnold SV et al. Health status after transcatheter mitral valve repair in heart failure and secondary mitral regurgitation. *JACC* Mar 2019, 25951; DOI: 10.1016/j.jacc.2019.02.010. 14. Obadia JF et al. Percutaneous repair or medical treatment for secondary mitral regurgitation. *N Engl J Med* 2018; 379:2297-2306. 15. Stone G. Pivotal transcatheter FMR device trials: Focus on COAPT and MITRA-FR, with implications for other transcatheter mitral valve device investigations. FDA Townhall, TCT 2018. 16. Glower DD, Kar S, Trento A, Lim DS, Bajwa T, Quesada R, et al. Percutaneous mitral valve repair for mitral regurgitation in high risk patients. *JACC*. 2014;64(2):172-81. 17. Whitlow PL, Feldman T, Pedersen WR, Lim DS, Kipperman R, Smalling R, et al. Acute and 12-month results with catheter based mitral valve leaflet repair: the EVEREST II (Endovascular Valve Edge-to-Edge Repair) High Risk Study. *JACC*. 2012 Jan 10;59(2):130-9. PubMed PMID: 22222076. 18. Maisano F, Franzen O, Baldus S, et al. Percutaneous mitral valve interventions in the real world: early and 1-year results from the ACCESS-EU, a prospective, multicentre, nonrandomized post-approval study of the MitraClip Therapy in Europe. *JACC*. 2013;62(12):1052-61. 19. Butter C, Franzen O, Treede H, et al. Analysis of outcomes at 12 months by mitral regurgitation etiology in the ACCESS-EUROPE phase I study. In: EuroIntervention, editor. EuroPCR; Paris: EuroIntervention; 2013. 20. Hoffman I. Analysis of outcomes at 12 months by gender in the ACCESS-EUROPE Phase I Study. Presented at PCR London Valves 2013, London, UK. 2013. 21. Kar S, Lim D, Smalling R, Whisenant B, Rammohan C, Fail PS, et al. The EVEREST II REALISM continued access study: effectiveness of transcatheter reduction of significant mitral regurgitation in surgical candidates. *JACC*. 2013;61(10):E1959. 22. Lim DS, Reynolds MR, Feldman T, Kar S, Herrmann HC, Wang A, et al. Improved functional status and quality of life in prohibitive surgical risk patients with degenerative mitral regurgitation following transcatheter mitral valve repair with the MitraClip System. *JACC*. 2013 Oct 24. PubMed PMID: 24184254. 23. Puls M, Lubos E, Boekstegers P, et al. One-year outcomes and predictors of mortality after MitraClip therapy in contemporary clinical practice: results from the German transcatheter mitral valve interventions registry. *European Heart Journal*. 2016;37(8):703-12. 24. Wiebe J, Franke J, Lubos E, Boekstegers P, Schillinger W, Quarrak T, et al. Percutaneous mitral valve repair with the MitraClip system according to the predicted risk by the logistic EuroSCORE: Preliminary results from the German Transcatheter Mitral Valve Interventions (TRAMI) registry. *Catheterization and Cardiovascular Interventions*. 2014 Mar 24. PubMed PMID: 24664460. 25. Rudolph V, Hüntgeburth M, Von Bardeleben RS, et al. Clinical outcome of critically ill, not fully recompensated, patients undergoing MitraClip therapy. *European Journal Heart Failure*. 2014;16:1223-9.

**CAUTION:** This product is intended for use by or under the direction of a physician. Prior to use, reference the Instructions for Use provided inside the product carton (when available) or at [eifu.abbottvascular.com](http://eifu.abbottvascular.com) or at [medical.abbott/manuals](http://medical.abbott/manuals) for more detailed information on Indications, Contraindications, Warnings, Precautions and Adverse Events. Information contained herein for DISTRIBUTION in Australia and New Zealand ONLY. Check the regulatory status of the device in areas where CE marking is not the regulation in force.

Illustrations are artist's representations only and should not be considered as engineering drawings or photographs. Photos on file at Abbott.

Abbott Australasia Pty Ltd: 299 Lane Cove Road, Macquarie Park, NSW, Australia 2113, Tel: 1800 550 939

Abbott Laboratories (NZ) Limited: Building D, 4 Pacific Drive, Mount Wellington, Auckland 1641, Tel: 0800 827 285

™ indicates a trademark of the Abbott Group of Companies.

[www.Cardiovascular.Abbott](http://www.Cardiovascular.Abbott)

©2020 Abbott. All rights reserved. ANZ SH0143-EN 01/20 | Item approved for use in ANZ only.

