

# Gallant™ HF CRT-D

CDHFA500T



Compatible with myMerlinPulse™ app

## Product Highlights

- Bluetooth® Low Energy (LE) communication enabling smartphone connectivity through data encryption
- MultiPoint™ pacing delivers multiple LV pacing pulses per cardiac cycle in both LV only and BiV pacing modes
- SyncAV™ Plus CRT technology offers dynamic AV timing with adaptive programming to ensure BiV pacing with or without MultiPoint pacing
- Improved shape with reduced volume and thickness
- 40J delivered energy safety shock option for enhanced safety margin
- DeFT Response™ technology offers noninvasive programming options to optimize rescue therapy to each patient's unique physiology and changing conditions
- VF Therapy Assurance decreases time to treatment for arrhythmias in patients who are likely to be hemodynamically unstable
- Antitachycardia pacing (ATP) while charging and prior to charging in the VF zone further extends the programming options for terminating tachyarrhythmias without a high-voltage shock
- ShockGuard™ technology with DecisionTx™ programming designed to reduce inappropriate therapy and minimize the need for programming adjustments at implant
  - SecureSense™ RV lead noise discrimination algorithm detects sustained lead noise and short bursts of oversensing that would otherwise go unnoticed or potentially lead to one or more inappropriate shocks
  - Far Field MD™ morphology discrimination and chamber onset discrimination enhances SVT and VT discrimination for reduced inappropriate therapies
- SenseAbility™ sensing algorithm feature provides the flexibility to fine-tune programming around T-wave oversensing without decreasing sensitivity
- The Gallant HF CRT-D and Quartet™ quadripolar LV lead feature four pacing electrodes and 13 pacing vectors to provide more options and greater control to address implant complications such as diaphragmatic stimulation and high pacing thresholds
- Easily test and program with Auto VectSelect Quartet™ multivector testing, offering an efficient workflow for complete results and programming
- DynamicTx™ over-current detection algorithm automatically changes shock configurations to ensure delivery of high-voltage therapy when high current is detected
- MRI-Ready device tested in combination with MR Conditional leads for full-body scans using a 1.5T or 3T (Tesla) field strength MRI Scanner\*
- The CorVue™ thoracic impedance feature measures transthoracic impedance changes over time to provide additional insight into the patient's heart failure condition
- Cold can programmability provides an additional RV-SVC shock configuration to decouple the can from the shocking vector parameters in cases of lead problems
- Premature Atrial Contraction (PAC) Response to avoid pacing the atrium in a vulnerable zone
- Physiologic rate responsive AV Delay and PVARP
- QuickOpt™ timing cycle optimization provides quick and effective optimization at the push of a button
- Dual patient notification: audio notification through the device and visual notification via myMerlinPulse app

## Ordering Information

Contents: Cardiac Pulse Generator

| MODEL NUMBER | DIMENSIONS (L × W × H) (MM) | WEIGHT (G) | VOLUME (CC) | CONNECTOR DEFIBRILLATION | CONNECTOR SENSE/PACE | CONNECTOR SENSE/PACE - LEFT VENTRICLE |
|--------------|-----------------------------|------------|-------------|--------------------------|----------------------|---------------------------------------|
| CDHFA500T    | 79 × 51 × 12                | 81         | 37          | DF-1                     | IS-1                 | IS-4                                  |

\*See MRI Scan Parameters in MRI-Ready Systems Manual.

## Product Specifications

| PARAMETER SPECIFICATIONS                                      |   |
|---|---|
| Model   | CDHFA500T   |
| Telemetry   | Bluetooth® LE Communication   |
| Delivered/Stored Energy                                       | 40/45 J   |
| Volume  | 37 cc   |
| Weight  | 81 g  |
| Size  | 79 × 51 × 12 mm   |
| Defibrillation Lead Connection                                | DF-1  |
| Atrial Sense/<br>Pace Lead Connection                         | IS-1 in-line bipolar  |
| Ventricular Sense/<br>Pace Lead Connection                    | IS-1 in-line bipolar  |
| Left Ventricular<br>Pace Lead Connection                      | IS4-LLLL  |
| High Voltage Can  | Electrically active titanium can  |
| Parameter   | Settings  |
| <b>Biventricular Pacing</b>                                   |   |
| VectSelect Quartet™<br>Programmable<br>LV Pulse Configuration | Distal Tip 1 - Mid 2; Distal Tip 1 - Proximal 4; Distal Tip 1 - Mid 3; Distal Tip 1 - RV Coil; Mid 2 - Mid 3; Mid 2 - Proximal 4; Mid 2 - RV Coil; Mid 3 - Mid 2; Mid 3 - Proximal 4; Mid 3 - RV Coil; Proximal 4 - Mid 2; Proximal 4 - Mid 3; Proximal 4 - RV Coil |
| MultiPoint™ Pacing  | LV1, LV2  |
| Delay MultiPoint Pacing                                       | Delay 1: 5; 10; ... 80 ms<br>Delay 2: 5; 10; ... 50 ms  |
| V. Triggering   | On; Off   |
| QuickOpt™ Timing Cycle<br>Optimization                        | Sensed/Paced AV delay, Interventricular pace delay  |
| V-V Timing  | Simultaneous <sup>1</sup> ; RV First; LV First  |
| Interventricular Pace Delay                                   | RV First 10–80/LV First 15–80 ms  |
| Ventricular Sensing   | RV only (not programmable)  |
| Ventricular Pacing Chamber                                    | RV only; LV only; Biventricular   |
| SyncAV™ Plus CRT<br>Technology Delta                          | If Type = Percentage: -10; - 15;...-70%<br>If Type = Fixed: -10; -20;...-120 ms; Off  |
| MPP PVAB  | 125-260 ms  |
| <b>AF Management</b>  |   |
| AF Suppression™ Pacing  | On; Off   |
| No. of Overdrive Pacing Cycles                                | 15-40   |
| Maximum AF Suppression Rate                                   | 80-150 bpm  |

## Product Specifications

| PARAMETER SPECIFICATIONS                            |   |
|---|---|
| <b>Sensing/Detection</b>                            |   |
| SenseAbility™ Sensing Algorithm                     | Automatic sensitivity control adjustment for atrial and ventricular events  |
| Low Frequency Attenuation                           | On; Off   |
| Threshold Start                                     | Post-Sensed: 50; 62.5; 75; 100%<br>Post-Paced, Atrial: 0.2-3.0 mV<br>Post-Paced, Ventricular: Auto, 0.2-3.0 mV  |
| Decay Delay   | Post-Sensed: 0-220 ms<br>Post-Paced, Atrial: 0-220 ms<br>Post-Paced, Ventricular: Auto, 0-220 ms  |
| Ventricular Sense Refractory                        | 125; 157 ms   |
| Detection Zones                                     | 3 zone programming – 1 zone; 2 zones; or 3 zones (VT-1; VT-2; VF)   |
| SVT Discriminators                                  | AV Rate Branch; Arrhythmia Onset (Chamber Onset or Sudden Onset); Interval Stability; AV Association Morphology; Discrimination (Far Field MD™ Morphology Discrimination or Original MD) with Automatic Template Update                   |
| Monitor Mode  | Detection, discrimination, and diagnostics, no therapy delivery (VT or VT-1 zone)   |
| Discrimination Modes                                | On; Passive; Off  |
| SVT Upper Limit                                     | 150-240 bpm   |
| SVT Discrimination Timeout                          | 20s-60 min; Off   |
| Reconfirmation                                      | Continuous sensing during charging  |
| SecureSense™ RV Lead Noise Discrimination Algorithm | On; On with Timeout; Passive; Off   |
| VF Therapy Assurance                                | On; Off   |
| <b>Antitachycardia Pacing Therapy</b>               |   |
| ATP Configurations                                  | Ramp; Burst; Scan; 1 or 2 schemes per VT zone   |
| ATP in VF Zone                                      | ATP While Charging; ATP Prior to Charging; Off  |
| ATP Upper Rate Cutoff                               | 150–300 bpm   |
| Burst Cycle Length                                  | Adaptive (50%-100%); Fixed (200-550 ms)   |
| Min. Burst Cycle Length                             | 150-400 in increments of 5 ms   |
| Readaptive  | On; Off   |
| Number of Bursts/Stimuli                            | 1-15 with 2–20 Stimuli  |
| Add Stimuli per Burst                               | On; Off   |
| ATP Pulse Amplitude                                 | 7.5 V independent from bradycardia and post-therapy pacing  |
| ATP Pulse Width                                     | 1.0 or 1.5 ms independently programmable from bradycardia and post-therapy pacing   |
| <b>High-Voltage Therapy</b>                         |   |
| DynamicTx™ Over-current Detection Algorithm         | On; Off   |
| DeFT Response™ Technology                           | Programmable pulse width for P1/P2 and tilt   |
| High-Voltage Output Mode                            | Fixed Pulse Width; Fixed Tilt   |
| Waveform  | Biphasic; Monophasic  |
| RV Polarity   | Cathode (-); Anode (+)  |
| Electrode Configuration                             | RV to Can; RV to SVC/Can; RV to SVC   |
| <b>Bradycardia Pacing</b>                           |   |
| Permanent Modes                                     | DDD(R); DDT(R); DDI(R); VVT(R); VVI(R); AAI(R); Off   |
| Temporary Modes                                     | DDD; DDT; DDI; VVT; VVI; AAI; AAT; DOO; VOO; AOO; Off   |
| Activity Sensor                                     | On; Passive; Off  |
| Programmable Rate and Delay Parameters              | Base Rate (bpm); Rest Rate (bpm); Maximum Tracking Rate (bpm); Max Trigger Rate (bpm); Maximum Sensor Rate (bpm); Paced AV Delay (ms); Sensed AV Delay (ms); Rate Responsive AV Delay; Hysteresis Rate (bpm); Rate Hysteresis with Search |
| Pulse Amplitude                                     | 0.25-7.5 V  |
| Pulse Width   | 0.05, 0.1-1.5 ms  |
| LVCap™ Confirm Feature, LVCap™ 2 Confirm Feature    | Setup; On; Monitor; Off   |
| RVCap™ Confirm Feature                              | Setup; On; Monitor; Off   |

## Product Specifications

| PARAMETER SPECIFICATIONS   |   |
|--|---|
| <b>Bradycardia Pacing</b>  |   |
| ACap™ Confirm Feature  | On; Monitor; Off  |
| Auto Mode Switch (AMS)   | DDI(R); DDT(R); VVI(R); VVT(R); Off   |
| Atrial Tachycardia Detection Rate  | 110-300 bpm   |
| AMS Base Rate  | 40; 45; ... 135 bpm   |
| PMT Detection/Termination  | Atrial Pace; Passive; Off   |
| Rate Responsive PVARP  | Low; Medium; High; Off  |
| Rate Responsive V Pace Refractory  | On; Off   |
| PAC Response   | On; Off   |
| PAC Response interval  | 200-400 ms  |
| Shortest AV Delay  | 25-120 ms   |
| <b>Post-Therapy Pacing (Independently programmable from Bradycardia and ATP)</b> |   |
| Post-Shock Pacing Mode   | AAI; VVI; DDI; or DDD; Off  |
| Post-Shock Base Rate   | 30-100 bpm  |
| Post-Shock Pacing Duration   | 0.5; 1; 2.5; 5; 7.5; or 10 min; Off   |
| <b>Device Testing/Induction Methods</b>  |   |
| DC Fibber™ Induction Method Pulse Duration                                       | 0.5-5.0 sec   |
| Burst Fibber Cycle Length  | 20-100 ms   |
| Noninvasive Programmed Stimulation (NIPS)  | 2-25 stimuli with up to 3 extra stimuli   |
| <b>Patient Notifiers</b>   |   |
| Programmable Notifiers (On; Off)   | BatteryAssurance™ alert; Possible HV circuit damage; HV charge timeout; Long charge time for Capacitor Maintenance; Device at ERI; Atrial pacing lead impedance out of range. Ventricular pacing lead impedance out of range; High-voltage lead impedance out of range; AT/AF Episode duration; AT/AF Burden; High ventricular rate during AT/AF; SecureSense™ lead noise detection; Non-sustained ventricular oversensing; Ventricular pacing percentage greater than limit. |
| Device Parameter Reset   | On  |
| Entry into Backup VVI Mode   | On  |
| Auditory Duration  | 2; 4; 6; 8; 10; 12; 14; 16 sec  |
| Number of Audio alerts per Notification  | 2   |
| Number of Notifications  | 1-16  |
| Time Between Notifications   | 10; 22 hours  |
| <b>Electrograms and Diagnostics</b>  |   |
| Stored Electrograms  | 30 minutes (2 user programmable + discrimination channel), up to 1 minute programmable pre-trigger data per VT/VF electrograms; additional triggers include lead noise detection, non-sustained ventricular oversensing, morphology template updates, atrial episode, PMT termination, PAC response, magnet reversion, noise reversion  |
| Therapy Summary  | Diagram of therapies delivered  |
| Episodes Summary   | Directory listing of up to 60 episodes with access to more details including stored electrograms  |
| Lifetime Diagnostics   | History of bradycardia events and device-initiated charging   |
| AT/AF Burden Trend   | Trend data and counts   |
| Ventricular HV Lead Impedance  | Multi-Vector Trend Data   |
| Histograms and Trends  | Event Histogram; AV Interval Histogram; Mode Switch or AT/AF Duration Histogram; Peak Filtered Atrial Rate during Atrial Arrhythmia Histogram; Atrial Heart Rate Histogram; Ventricular Heart Rate Histogram; AT/AF Burden; Exercise and Activity Trending; V Rates during AMS; DirectTrend™ reports up to 1 year   |
| PMT Data   | Information regarding PMT detections  |
| Real-Time Measurements (RTM)   | Pacing lead impedances; High-voltage lead impedances; Signal amplitudes   |
| CorVue Thoracic Impedance  | On; Off   |
| CorVue Thoracic Impedance Threshold  | 8-18 days   |

## Product Specifications

| PARAMETER SPECIFICATIONS          |   |
|-----------------------------------|---|
| MRI Settings                      |   |
| Tachy Therapy                     | Disabled  |
| MRI Mode                          | DOO; VOO; AOO; Pacing Off                                     |
| MRI Base Rate                     | 30 - 100 bpm  |
| MRI Paced AV Delay                | 25-110 ms   |
| MRI RA and RV Pulse Amplitude     | 5.0 or 7.5 V  |
| MRI RA and RV Pulse Width         | 1.0 ms  |
| MRI RA and RV Pulse Configuration | Bipolar   |
| MRI LV Pulse Configuration        | D1-M2, D1-M3, D1-P4, M2-M3, M2-P4, M3-M2, M3-P4, P4-M2, P4-M3 |
| MRI LV Pulse Amplitude            | 0.25-7.5 V  |
| MRI LV Pulse Width                | 0.05-1.5 ms   |
| MRI V Pacing Chamber              | RV Only, LV+RV (Simultaneous)                                 |
| MRI Timeout                       | 3; 6; 9; 12; 24 hours; Off                                    |

| MRI SCAN PARAMETERS <sup>§</sup>   |                |                             |             |
|------------------------------------|----------------|-----------------------------|-------------|
| Lead Model                         | Magnet (Tesla) | RF Transmit Conditions      | Scan Region |
| <b>Durata™ Defibrillation Lead</b> | 1.5 T / 3 T    | Normal<br>Operating<br>Mode | Full-body   |
| 7120 (lead lengths: 65 cm)         |                |                             |             |
| 7122 (lead lengths: 60, 65 cm)     |                |                             |             |
| <b>Optisure™ Lead</b>              | 1.5 T / 3 T    |                             |             |
| LDA210 (lead lengths: 65 cm)       |                |                             |             |
| <b>Quartet™ LV Lead</b>            | 1.5 T / 3 T    |                             |             |
| 1456Q (lead length: 86 cm)         |                |                             |             |
| 1457Q (lead length: 86 cm)         |                |                             |             |
| 1458Q (lead length: 86 cm)         |                |                             |             |
| 1458QL (lead length: 86 cm)        | 1.5 T / 3 T    |                             |             |
| <b>Tendril™ STS Pacing Lead</b>    |                |                             |             |
| 2088TC (lead lengths: 46, 52 cm)   | 1.5 T / 3 T    |                             |             |
| <b>UltiPace™ Pacing Lead</b>       |                |                             |             |
| LPA1231 (Lead lengths: 46, 52 cm)  | 1.5 T / 3 T    |                             |             |

<sup>†</sup> LV first with 10 ms interventricular delay.

<sup>§</sup> For additional information about specific MR Conditional CRT-Ds and leads, including scan parameters, warnings, precautions, adverse conditions to MRI scanning, and potential adverse events, please refer to the Abbott MRI-Ready Systems Manual at [manuals.eifu.abbott.com](https://manuals.eifu.abbott.com).

**Brief Summary:** This product is intended for use by or under the direction of a Physician. Prior to using these devices, please review the Instructions for Use for a complete listing of indications, contraindications, warnings, precautions, potential adverse events and directions for use.

**Abbott**  
15900 Valley View Court, Sylmar, CA 91342  
Tel: +1 818 362 6822  
Abbott.com

™ Indicates a trademark of the Abbott group of companies.

‡ Indicates a third party trademark, which is property of its respective owner.

© 2024 Abbott. All Rights Reserved.  
MAT-2407279 v1.0 | Item approved for OUS use.

