

DETECTING ATRIAL FIBRILLATION

WITH ASSERT-IQ™ INSERTABLE CARDIAC MONITOR (ICM)



ATRIAL FIBRILLATION (AFIB) STATISTICS

AFIB IS **THE MOST COMMON TYPE** OF CARDIAC ARRHYTHMIA.¹



33.5 MILLION PEOPLE
worldwide have AFib, which is
1.5-2% of the world's population.²

AFIB RELATION TO STROKE



Patients with AFib are **5x MORE LIKELY**
to have a stroke.^{1,2}

**30-
40%**

Approximately **30-40%** of ischemic stroke patients are diagnosed with cryptogenic (unexplained) stroke following diagnostic workup. Commonly referred to as CS.³⁻⁵



AFIB DETECTION WITHIN CRYPTOGENIC STROKE PATIENTS



Long-term monitoring reveals AFib in **~30%** of cryptogenic stroke (CS) patients.⁶

LONG-TERM MONITORING WITH AN INSERTABLE CARDIAC MONITOR (ICM)

**3X
FASTER**

ICMs **DETECT AFIB 3X FASTER** than external cardiac monitoring methods for CS patients.⁷



Actual size represented




ICMs DETECT MORE AFIB THAN CONVENTIONAL METHODS

Cryptogenic Stroke and Underlying Atrial Fibrillation Study (CRYSTAL AF)⁸

The CRYSTAL AF study evaluated the incidence of paroxysmal AFib in patients with cryptogenic stroke and transient ischemic attack (TIA) using an ICM.

Detection of AFib by 36 Months



PATIENT IMPACT

84 DAYS

is the median time to AFib detection in cryptogenic stroke patients.⁸


HIGHER DIAGNOSIS RATES

Longer monitoring duration leads to higher AFib diagnosis rates.⁸

FASTER AFIB DIAGNOSIS FOUND WITH REAL-WORLD EVIDENCE

A nationwide study of U.S. Medicare beneficiaries with cryptogenic stroke who received an ICM or an external cardiac monitor (ECM) after ischemic stroke demonstrates ICMs provide effective AFib detection, with almost three times more AFib detected compared to ECMs.⁷

New AFib Diagnosis in ICM vs ECM Patients



3X FASTER

AFib diagnosis for ischemic stroke patients monitored with an Abbott ICM compared to ECMs.⁷

~3X MORE

ICM patients prescribed with an OAC at 1.5 years compared to ECM patients.⁷

GUIDELINES RECOMMEND AN ICM FOR AFIB DETECTION FOR PATIENTS WITH CRYPTOGENIC STROKE

2019 AHA/ACC/HRS Atrial Fibrillation Guidelines

Recommends use of an ICM in patients with cryptogenic stroke (Class IIa, LOE B-R).⁹

RECOMMENDATION	COR	LOE
In patients with cryptogenic stroke (i.e., stroke of unknown cause) in whom external ambulatory monitoring is inconclusive, implantation of a cardiac monitor (loop recorder) is reasonable to optimize detection of silent atrial fibrillation.	IIa	B-R

2020 European Society Guidelines: Management of Atrial Fibrillation

In patients with stroke and without diagnosed AFib, additional monitoring by non-invasive ECG monitor or injectable recorder should be considered in the long-term.¹⁰

RECOMMENDATION	CLASS	LEVEL
In stroke patients, additional ECG monitoring by long-term non-invasive ECG monitors or implanted loop recorders should be considered to document silent atrial fibrillation.	IIa	B

COST EFFECTIVENESS OF CONTINUOUS MONITORING

- Cost-effectiveness models demonstrate ICMs implanted immediately versus delayed (after Holter monitoring) save payers money per quality-adjusted life year (QALY).
 - **\$26,342 to \$42,967 savings** per QALY with an ICM.¹¹⁻¹³
 - Cost savings increase among subgroup analyses with lower CHA2DS2-VASc scores.^{11,12}
- In the U.S. model, immediate ICM implant achieves higher rates of AFib detection and treatment, compared to standard of care intermittent ECG.¹¹



ADVANCED AFIB ALGORITHM

- With Assert-IQ™ ICM, you can detect arrhythmias more accurately,* with industry-leading advanced algorithms that provide clinically actionable data.
- **5-step AFib detection discriminator** focuses on R-R interval patterns and P-waves in EGMs to verify if a detected AFib event is true or false.

5-STEP SEQUENTIAL AFIB DETECTION DISCRIMINATOR

Interval-Based Steps

1. Pattern Recognition
2. R-Wave Undersensing
3. T-Wave Oversensing

Morphology-Based Steps

4. P-Wave Oversensing
5. P-Wave Detection



**Compared to predicate ICM devices*

WHY ASSERT-IQ™ ICM?

SMALL

Slim design; not noticeable after insertion in most patients.

INTELLIGENT

Advanced diagnostic capabilities enable more informed decision making.*

FAST

The transmitter checks in with the device every 20 seconds, so patients are continuously monitored.

LATEST BLUETOOTH®-ENABLED TECHNOLOGY

Automatic, seamless transmission of patient information while they go about their normal day.

LONGEST BATTERY LIFE

Building on a reputation for innovation, Assert-IQ EL+ is the longest lasting Bluetooth®-enabled ICM that operates at full functionality beyond six years with no compromises,¹⁴⁻²² making it an effective tool to support the expanding patient management needs of your practice.**

* Compared to predicate ICM devices

** DM5500

NO COMPROMISE WITH **MRI SCANS**

Assert-IQ™ ICM is MR Conditional for both 1.5 Tesla (T) and 3T full-body scans. There is no wait limit after insertion, which allows timely patient care.*



Choose between two longevity options to support atrial fibrillation detection and arrhythmia management for your patients.

Assert-IQ™ 3/3+

(DM5000/DM5300)

3 YEAR; May be preferred for traditional cardiac monitoring.



Assert-IQ™ EL+

(DM5500)

6 YEAR; Longer term monitoring; May be preferred for AFib management.



Actual sizes represented

Assert-IQ™ ICM. Only from Abbott.

Powered by SyncUP™ Remote Monitoring Support and myMerlin™ Mobile App, Assert-IQ ICM provides the flexibility to deliver data the way you and your patients need.

* No minimum length of time requirement between implant and MRI scan. For additional information, including warnings, precautions, and potential adverse events, please refer to the Abbott MRI-Ready Systems Manual at <https://medical.abbott/manuals> or visit cardiovascular.abbott/mrireaddy.

References

- Centers for Disease Control and Prevention. About Heart Disease - Atrial Fibrillation. Updated October 14, 2022. Accessed March 28, 2024. <https://www.cdc.gov/heartdisease/atrial-fibrillation.htm>
- Chugh SS, Havmoeller R, Narayanan K, et al. Worldwide Epidemiology of Atrial Fibrillation: A Global Burden of Disease 2010 Study. *Circulation*. 2014;129(8):837-847. doi:10.1161/CIRCULATIONAHA.113.005119.
- Grau AJ, Weimar C, Buggle F, et al. Risk Factors, Outcome, and Treatment in Subtypes of Ischemic Stroke: The German Stroke Data Bank. *Stroke* 2001; 32: 2559-2566.
- Adams HP Jr, Bendixen BH, Kappelle LJ, et al. Classification of Subtype of Acute Ischemic Stroke. Definitions for Use in a Multicenter Clinical Trial. TOAST. Trial of Org 10172 in Acute Stroke Treatment. *Stroke* 1993; 24: 35-41.
- Kolominsky-Rabas PL, Weber M, Gefeller O, Neundorfer B, Heuschmann PU. Epidemiology of Ischemic Stroke Subtypes According to TOAST Criteria: Incidence, Recurrence, and Long-Term Survival in Ischemic Stroke Subtypes: A Population-Based Study. *Stroke* 2001;32:2735.
- Kloosterman EM, Rosman JJ, Berkowitz EJ, Rosenbaum M, Wittenstein ZA. Treatment of Cryptogenic Stroke Patients with Atrial Fibrillation Detected by Insertable Cardiac Monitors Reduces Recurrent Stroke Risk to Background Levels. *The Journal of Innovation in Cardiac Rhythm Management*. 2021;12(12):4812-4817. doi:10.19102/icrm.2021.121204
- Boriani G, Auricchio A, Botto GL, et al. Insertable Cardiac Monitoring Results in Higher Rates of Atrial Fibrillation Diagnosis and Oral Anticoagulation Prescription after Ischaemic Stroke. *EP Europace*. 2023; eua212, <https://doi.org/10.1093/europace/euad212>
- Sanna T, Diener H-C, Passman RS, et al. Cryptogenic Stroke and Underlying Atrial Fibrillation. *The New England Journal of Medicine*. 2014;370(26):2478-2486. doi:10.1056/NEJMoa1313600
- January CT, Wann LS, Calkins H, et al. 2019 AHA/ ACC/HRS Focused Update of the 2014 AHA/ACC/ HRS Guideline for the Management of Patients with Atrial Fibrillation: A Report of the American College of Cardiology/American Heart Association Task Force on Clinical Practice Guidelines and the Heart Rhythm Society. *Heart Rhythm*. 2019;16(8):e66-e93. doi:10.1016/j.hrthm.2019.01.024.
- G. Hindricks, T. Potpara, N. Dagres, et al. 2020 ESC Guidelines for the Diagnosis and Management of Atrial Fibrillation Developed in Collaboration with the European Association for Cardio-Thoracic Surgery (EACTS). 2021;26(9). doi:10.15829/1560-4071-2021-4701.
- Sawyer LM, Witte KK, Reynolds MR, et al. Cost-Effectiveness of an Insertable Cardiac Monitor to Detect Atrial Fibrillation in Patients with Cryptogenic Stroke. *Journal of Comparative Effectiveness Research*. 2021;10(2):127-141. doi:10.2217/cer-2020-0224.
- Edwards SJ, Wakefield V, Jhita T, Kew K, Cain P, Marcenik G. Implantable Cardiac Monitors to Detect Atrial Fibrillation after Cryptogenic Stroke: A Systematic Review and Economic Evaluation. *Health Technology Assessment (Winchester, England)*. 2020;24(5):1-184. doi:10.3310/hta24050.
- Vincent Thijs, Klaus K. Witte, Carmel Guarnieri, et al. Cost-effectiveness of Insertable Cardiac Monitors for Diagnosis of Atrial Fibrillation in Cryptogenic Stroke in Australia. *Journal of Arrhythmia*. 2021;37(4):1077-1085. doi:10.1002/joa3.12586
- Abbott. Assert-IQ™ ICM User Manual.
- Medtronic. REVEAL LINQ™ LNQi1 Insertable Cardiac Monitor and Patient Assistant PA96000 Clinician Manual. Updated August 26, 2015. Accessed February 23, 2024. https://www.medtronic.com/content/dam/emanuals/crdm/CONTRIB_215651.pdf
- Medtronic. LINQ II™ LNQi22 Insertable Cardiac Monitor Clinician Manual. Updated September 01, 2022. Accessed February 23, 2024. https://www.medtronic.com/content/dam/emanuals/crdm/M032283C001B_view.pdf
- Boston Scientific. User's Manual, LUX-Dx⁺ Insertable Cardiac Monitor System M301, 2925, 2935. Updated April 2023. Accessed February 23, 2024. https://www.bostonscientific.com/content/dam/elabeling/crm/51656369-001_LUX-Dx_ICM_UM_en_S.pdf
- Boston Scientific. User's Manual, LUX-Dx II/II+⁺ Insertable Cardiac Monitor System M302, M312, 2925, 2929, 2935, 2939. Updated August 2023. Accessed February 23, 2024. https://www.bostonscientific.com/content/dam/elabeling/crm/51583079-001_LUX-Dx_ICM_UM_en_S.pdf
- Biotronik. Technical Manual BIOMONITOR III⁺. Updated December 10, 2020. Accessed February 23, 2024. https://manuals.biotronik.com/emanuals-professionals/?country=US&product=ImplCardMon/BioMonitor3/BioMonitor3_US
- Biotronik. Technical Manual BIOMONITOR IIIIm⁺. Updated December 10, 2020. Accessed February 23, 2024. https://manuals.biotronik.com/emanuals-professionals/?country=US&product=ImplCardMon/BioMonitor3m/BioMonitor3m_US
- Biotronik. Technical Manual BIOMONITOR IV⁺. Updated July 17, 2023. Accessed February 23, 2024. https://manuals.biotronik.com/emanuals-professionals/?country=US&product=ImplCardMon/BioMonitor4/BioMonitor4_US

Data on File. Abbott - Report 90984075.

Rx Only

Brief Summary: Prior to using these devices, please review the Instructions for Use for a complete listing of indications, contraindications, warnings, precautions, potential adverse events and directions for use.

Indications for Use: The Assert-IQ™ ICM is indicated for the monitoring and diagnostic evaluation of patients who experience unexplained symptoms that may be cardiac-related such as: dizziness, palpitations, chest pain, syncope, and shortness of breath, as well as patients who are at risk for cardiac arrhythmias such as bradycardia, tachycardia, and sinus pauses. The Assert-IQ ICM is also indicated for patients who have been previously diagnosed with atrial fibrillation (AF) or who are susceptible to developing AF. The Assert-IQ ICM is intended to be inserted subcutaneously in the left pectoral region, also described as the left anterior chest wall. The Assert-IQ ICM has not been specifically tested for pediatric use.

Intended Use: The Assert-IQ ICM is intended to help physicians and clinicians monitor, diagnose and document the heart rhythm in patients who are susceptible to cardiac arrhythmias and unexplained symptoms by detecting arrhythmias and transmitting data for review.

Contraindications: There are no known contraindications for the insertion of the Assert-IQ ICM. However, the patient's particular medical condition may dictate whether or not a subcutaneous, chronically inserted device can be tolerated.

Potential Adverse Events: Possible adverse events (in alphabetical order) associated with the device, include the following: allergic reaction, bleeding, chronic nerve damage, erosion, excessive fibrotic tissue growth, extrusion, formation of hematomas or cysts, infection, keloid formation and migration.

Refer to the User's Manual for detailed indications for use, contraindications, warnings, precautions and potential adverse events.

An Abbott mobile transmitter is available for patients without their own compatible mobile device.

™ Indicates a trademark of the Abbott group of companies.

† Indicates a third-party trademark, which is property of its respective owner.

Bluetooth and Bluetooth logo are registered trademarks of Bluetooth SIG, Inc.

Abbott

15900 Valley View Court, Sylmar, CA 91342
Tel: +1 818 362 6822
Abbott.com

© 2024 Abbott. All Rights Reserved.
MAT-2200807 v4.0 | Item approved for U.S. use.



Experience a Faster Diagnosis
with **Assert-IQ™ ICM**

