Assert-IQ™

Insertable Cardiac Monitor (ICM)

# Assert-IQ™ ICM

## Clinical Compendium

# CONTENTS

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## 数字で見るAssert-IQ™ ICM

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
Enhanced Algorithms Maintain  
AF Time-to-Diagnosis

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## ADVANCED ALGORITHMS

# Evaluating the impact of New Arrhythmia Detection Algorithms in an ICM

SHEHATA ET AL. APHRS POSTER 2022<sup>4</sup> (CONTINUED)

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


Abbott

# Evaluation of A Novel Premature Ventricular Contraction (PVC) Detection Algorithm in an ICM

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


|                           | Training | Testing |
|---------------------------|----------|---------|
| Sensitivity               | 78.6%    | 75.8%   |
| Specificity               | 99.5%    | 99.7%   |
| Positive Predictive Value | 70.8%    | 82.7%   |
| Negative Predictive Value | 99.7%    | 99.5%   |

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
## Evaluation of A Novel Premature Ventricular Contraction (PVC) Detection Algorithm in an ICM

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# ICM P-Wave Visibility in a New Clinical Report

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
| Post-implant time | Total beats analyzed (n) | Beats with visible P-wave (%) | EGMs with visible P-wave in >50% of beats (%) |
|-------------------|--------------------------|-------------------------------|---|
| Day 30            | 1429                     | 89.1%                         | 90.4%   |
| Day 60            | 1459                     | 90.8%                         | 92.0%   |
| Day 90            | 1500                     | 89.4%                         | 90.3%   |

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# REFERENCES

1. Gopinathannair R, Lakkireddy D, Manyam H, et al. Improving the Specificity of Atrial Fibrillation and Tachycardia Detection in an Insertable Cardiac Monitor. Presented at: Heart Rhythm Society (HRS); April 29 - May 1, 2022; San Francisco, CA.
2. Afzal MR, Gopinathannair R, Manyam H, et al. Development and Evaluation of A New Algorithm Enhancement to Improve Specificity of Pause Detection in An Insertable Cardiac Monitor. Presented at: Heart Rhythm Society (HRS); April 29 - May 1, 2022; San Francisco, CA.
3. Gardner RS, Quartieri F, Betts TR, et al. Reducing the Electrogram Review Burden Imposed by Insertable Cardiac Monitors. *J Cardiovascular Electrophysiology*. 2022;33(4):741-750. doi:10.1111/jce.15397.
4. Shehata MM, Manyam H, Gopinathannair R, et al. Evaluating the Impact of New Arrhythmia Detection Algorithms in an Insertable Cardiac Monitor. Presented at: Asia Pacific Heart Rhythm Society (APHRs); November 18-20, 2022; Singapore.
5. Data on File. Abbott - Report SJM-CFM-0919-0163.
6. Manyam H, Afzal MR, Gopinathannair R, et al. Evaluation of a Novel Premature Ventricular Contraction Detection Algorithm in an Insertable Cardiac Monitor. Presented at: Heart Rhythm Society (HRS); April 29 - May 1, 2022; San Francisco, CA.
7. Arnar DO, Mairesse GH, Boriani G, et al. Management of Asymptomatic Arrhythmias: a European Heart Rhythm Association (EHRA) Consensus Document, Endorsed by the Heart Failure Association (HFA), Heart Rhythm Society (HRS), Asia Pacific Heart Rhythm Society (APHRs), Cardiac Arrhythmia Society of Southern Africa (CASSA), and Latin America Heart Rhythm Society (LAHRS). *Europace : European Pacing, Arrhythmias, and Cardiac Electrophysiology : Journal of the Working Groups on Cardiac Pacing, Arrhythmias, and Cardiac Cellular Electrophysiology of the European Society of Cardiology*. 2019;21(6):844–845. doi:10.1093/europace/euz046
8. Shehata M, Nair DG, Qu F, et al. Insertable Cardiac Monitor P-wave Visibility in a New Clinical Report. Presented at: Asia Pacific Heart Rhythm Society (APHRs); November 18-20, 2022; Singapore.

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Bluetooth and Bluetooth logo are registered trademark of Bluetooth SIG, Inc. ™ Indicates a trademark of the Abbott group of companies.

©2024 Abbott. All rights reserved.

MAT-2403023 v1.0 | Item approved for Japan use only.

**Abbott**