GALLANT™ ICDs & CRT-Ds

Precision Shock Technology

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
Defibrillation Waveform Duration Adjustment Increases the Proportion of Acceptable Defibrillation Thresholds in Patients Implanted with Single-Coil Defibrillation Lead¹³

	A群：至適波形		B群：Tilt固定	
	NICM	ICM	NICM	ICM
14J	-	20 (59%)	-	24 (34%)
21J	7 (20%)	3 (9%)	17 (24%)	8 (11%)
25J	1 (3%)	3 (9%)	7 (10%)	3 (4%)
31J	0	0	6 (8%)	4 (6%)
35J	0	0	2 (3%)	0

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Benefit of millisecond waveform durations for patients with high defibrillation thresholds¹⁴

All (N=54)	11.0 ± 5.5J	8.8 ± 4.1J	P=0.0001
DFT ≥ 15J (N=13)	18.7 ± 4.1J	13.4 ± 3.5J	P=0.009

“Tuned” Defibrillation Waveforms Outperform 50/50% Tilt Defibrillation Waveforms: A Randomized Multi-Center Study¹⁵

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1. Improved ICD Shock Efficacy using Programmable Pulse Width. Michael Katcher, MD; Kevin Davis, BS; Nima Badie, PhD; Luke C. McSpadden, PhD; Ulrika Birgersdotter-Green, MD; Jacqueline Joza, MD; Robert Schaller, MD; Gabriel Mouchawar, PhD. Heart Rhythm Society Conference. April 24-27, 2025.
2. Blatt JA, Poole JE, Johnson GW, Callans DJ, Raitt MH, Reddy RK, Marchlinski FE, Yee R, Guarnieri T, Talajic M, Wilber DJ, Anderson J, Chung K, Wong WS, Mark DB, Lee KL, Bardy GH; SCD-HeFT Investigators. No benefit from defibrillation threshold testing in the SCD-HeFT (Sudden Cardiac Death in Heart Failure Trial). J Am Coll Cardiol. 2008 Aug 12;52(7):551-6.
3. Cha YM, Hayes DL, Asirvatham SJ, Powell BD, Cesario DA, Cao M, Gilliam FR 3rd, Jones PW, Jiang S, Saxon LA. Impact of shock energy and ventricular rhythm on the success of first shock therapy: the ALTITUDE first shock study. Heart Rhythm. 2013 May;10(5):702-8.
4. Wilkoff BL. et al., 2015 HRS/EHRA/APHRS/SOLAECE expert consensus statement on optimal implantable cardioverter-defibrillator programming and testing. Europace; 2016;18;159-183
5. Stiles MK. et al., 2019 HRS/EHRA/APHRS/LAHS focused update to 2015 expert consensus statement on optimal implantable cardioverter-defibrillator programming and testing. Heartrhythm 2019;17;1:e220-e228
6. Windecker S. et al., The Influence of Ventricular Fibrillation Duration on Defibrillation Efficacy Using Biphasic Waveforms in Humans. JACC 1999;33;1;33-8
7. Kroll MW et al., Fundam Clin Pharmacol 2010;24(5):561-573
8. Tokano et al. Long-Term Evaluation of the Ventricular Defibrillation Energy Requirement. J Cardiovasc Electrophysiol, Vol. 9, pp. 916-920. September 1998
9. Pelosi et al. Effect of Chronic Amiodarone Therapy on Defibrillation Energy Requirements in Humans. J Cardiovasc Electrophysiol. Vol. 11. pp. 736-740, July 2000
10. Fain ES. et al., Implantable Cardioverter Defibrillator: Ventritex® Cadence® JCE 1993;4;pp. 211-223
11. Kroll MW et al., Optimizing defibrillation waveforms for ICDs J Interv Card Electrophysiol(2007) 18;247-263
12. Kroll MW et al., A Minimal Model of the Single Capacitor Biphasic Defibrillation Waveform. PACE(1994); 17(Pt.I):1782-1792
13. Gabriels J. et al., Defibrillation Waveform Duration Adjustment Increases the Proportion of Acceptable Defibrillation Thresholds in Patients Implanted with Single-Coil Defibrillation Lead. Cardiology 2013;124:71-75
14. Denman RA et al., Benefit of millisecond waveform durations for patients with high defibrillation thresholds. Heart Rhythm 2006;3:536-541
15. Natarajan, S. et al., "Tuned" Defibrillation Waveforms Outperform 50/50° Tilt Defibrillation Waveforms: A Randomized Multi-Center Study. PACE 2007; 30:S139-S142

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TM Indicates a trademark of the Abbott group of companies.

©2026 Abbott. All rights reserved.

MAT-2601479 v1.0 | Item approved for Japan use only. 100179742