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Kathleen M. Chester

Sr. Vice President, Regulatory Affairs and Quality Assurance

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¹ Mehta D, Nayak H, Singson M, Chao S, Jorge E, Camunas L, Gomes J. Late Complications in Patients with Pectoral Defibrillator Implants with Transvenous Defibrillator Lead Systems: High Incidence of Insulation Breakdown. *PACE*, 1998; 21:1893-1900.

² De Lurgio D, Sathavorn C, Mera F, Leon A, Walter P, Langberg J. Incidence and Implications of Abrasion of Implantable Cardioverter-Defibrillator Leads. *American Journal of Cardiology*, 1997, 79:1409-1411.

³ Pavia S, Saliba W, Wilkoff B. Lead System Dysfunction, Diagnosis, and Therapy. In: Pacifico A. ed-in-chief. *Implantable Defibrillator Therapy: A Clinical Guide*. Norwell, Ma. Kluwer Academic Publishers 2002:259-278.

⁴ Kleemann T, Becker T, Doenges K, Vater M, Senges J, Schneider S, Saggau W, Weisse U, Seidl K. Annual Rate of Transvenous Defibrillation Lead Defects in Implantable Cardioverter-Defibrillators over a Period > 10 years. *Circulation*. 2007;115:2474-2480.

⁵ Blommaert D, Louagie Y, Eucher P, Collet B, Mancini I, Govarets G, Leclercq C, De Roy L. Are Cardioverter-Defibrillator Leads Reliable Over time? *NASPE Abstracts*, 2003, 26:1123.

⁶ Degeratu F, Khalighi K, Peters R, Shorofsky S, Gold M. Sensing Lead Failure in Implantable Defibrillators: A Comparison of Two Commonly Used Leads. *J Cardiovasc Electrophysio*, 2000, 11:21-24.

⁷ Valk S, Luijten R, Jordaens L. Insulation Damage in a Shock Wire: An Unexpected Fluoroscopic Image. *PACE*. 2010; 1-3.

⁸ Richards M, Warren C, Anderson M. Late failure of a single-coil Transvenous implantable cardioverter-defibrillator lead associated with conductor separation. *Europace Advanced Access* published March 27, 2010.

⁹ Jalal Z, Derval N, Ploux S, Bordachar P. Unusual failure of a multilumen, small-diameter implantable cardioverter-defibrillator lead. *Heart Rhythm*, 2009, Articles in Press.

¹⁰ Duray G, Israel C, Schmitt J, Hohnloser S. Implantable cardioverter-defibrillator lead disintegration at that level of the tricuspid valve. *Heart Rhythm Society*, 2008.

¹¹ HRS/EHRA Expert Consensus on Monitoring Cardiovascular Implantable Electronic Devices (CIED) dated April 2008

¹² Wilkoff et al., Transvenous Lead Extraction: Heart Rhythm Society Expert Consensus on Facilities, Training, Indications, and Patient Management, *Heart Rhythm*, Vol 6, No 7, July 2009.

¹³ Maisel et al., Recommendations from the Heart Rhythm Society Task Force on Lead Performance Policies and Guidelines, *Heart Rhythm*, Vol 6, No 6, June 2009.