FOCUS ON CLINICAL PERFORMANCE

Durata™ Lead Performance

Optim Lead Insulation

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

Optim® Lead Insulation Effects on SJM Tachycardia Lead Abrasion

Kaplan-Meier Analysis of US Returns Analysis Data

Year	1	2	3	4	5
Optim	147,866	100,578	60,435	28,900	9,159
Silicone	129,167	115,608	104,333	91,838	75,119

1 C. Jenney, J. Tan, A. Karicherla, J. Burke, and J. Helland, "A New Insulation Material for Cardiac Leads with Potential for Improved Performance," HRS 2005, HeartRhythm, 2, S318-S319 (2005).

2 J. Tan and C. Jenney, "Comparative In Vivo Biostability Study of A New Lead Insulation Material Versus Polyurethanes," HRS2006, Heart Rhythm, 3, S146 (2006).

3 T. Kleemann, T. Becker, K. Doenges, M. Vater, J. Senges, S. Schneider, W. Saggau, U. Weisse, and K. Seidl, "Annual Rate of Transvenous Defibrillation Lead Defects in Implantable Cardioverter-Defibrillators

Over a Period of >10 Years," Circulation, 115, 2474-2480 (2007).

ST. JUDE MEDICAL™

MORE CONTROL. LESS RISK.

St. Jude Medical is focused on reducing risk by continuously finding ways to put more control into the hands of those who save and enhance lives.

IMPLANTABLE ELECTRONIC SYSTEMS

CARDIOVASCULAR AND ABLATION TECHNOLOGIES

Global Headquarters

One St. Jude Medical Drive
St. Paul, Minnesota 55117
USA

+1 651 756 2000

+1 651 756 3301 Fax

**Implantable Electronic
Systems Division**

15900 Valley View Court
Sylmar, California 91342
USA

+1 818 362 6822

+1 818 364 5814 Fax

St. Jude Medical AB

Veddestavägen 19
175 84 Järfälla
Sweden

+46 8 474 40 00

+46 8 760 95 42 Fax

U.S. Division

6300 Bee Cave Road
Building Two, Suite 100
Austin, Texas 78746
USA

+1 512 732 7400

+1 512 732 2418 Fax

SJMprofessional.com

CAUTION: FEDERAL LAW (USA) RESTRICTS THIS DEVICE TO SALE, DISTRIBUTION AND USE BY OR ON THE ORDER OF A PHYSICIAN.
Brief Summary: Prior to using these devices, please review the User's Manual for a complete listing of indications, contraindications, warnings, precautions, potential adverse events and directions for use. Unless otherwise noted, ® or ™ indicates a registered or unregistered trademark or service mark owned by, or licensed to, St. Jude Medical, Inc. or one of its subsidiaries. ST. JUDE MEDICAL, the nine-squares symbol and MORE CONTROL. LESS RISK, are registered and unregistered trademarks and service marks of St. Jude Medical, Inc. and its related companies. ©2012 St. Jude Medical, Inc. All rights reserved.

ST. JUDE MEDICAL
MORE CONTROL. LESS RISK.