



myMerlin Mobile App Compatible Smart Phone and OS List

Popular mobile devices and operating systems are regularly tested to evaluate compatibility with the myMerlin mobile apps. We recommend checking this list before installing a new version on your phone or before using the app with a new phone. This table will be updated as other devices and operating systems (OS's) are evaluated. Jailbroken devices and beta OS versions are not supported.

Recommended app, devices, and OS versions:

App Version	Devices	Operating System
APP1001 V1.0.6000	iPhone 7, 8, 8+, X, XR, XS, SE, 11, 11 Pro, 11 Pro Max, 12	iOS 13 and greater
APP1000 V1.0.6000	Samsung Galaxy A01, A10e, A20e, A21s, A51, S8, S9, S10, S20, S21, S21Plus, J3, J7 Star	Android OS 9 and greater
	Google Pixel 2XL, 3, 3aXL, 4, 4XL, 5	
	Sony Xperia (802SO)	
	Sharp Aqueous Sense 3	
	LG Q70, V35 ThinQ, G6, K31	
	Motorola Moto G7, Z4, E6	
Xiaomi Mi 10T, Redmi Note 7		