



# AMPLATZER™ DUCT OCCLUDERS

## INDICATIONS

The Amplatzer™ Duct Occluder and Amplatzer™ Duct Occluder II are percutaneous, transcatheter occlusion devices intended for the nonsurgical closure of a patent ductus arteriosus (PDA) in patients with a weight of 6 kg or larger.


The Amplatzer Piccolo™ Occluder is a percutaneous, transcatheter occlusion device intended for the nonsurgical closure of a patent ductus arteriosus (PDA) in patients with a weight 700 grams and up at time of the procedure.



## MRI SAFETY INFORMATION

Through nonclinical testing, Amplatzer™ devices have been shown to be MR Conditional. Refer to the appropriate Instructions for Use to obtain more detailed MRI scanning information.

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the device. Therefore, optimization of MR imaging parameters to compensate for the presence of this device may be necessary.




### DEVICE SPECIFICATIONS AND DELIVERY SYSTEM COMPATIBILITY

#### Amplatzer™ Duct Occluder

#### Amplatzer™ TorqVue™ 180° Delivery System

Model/ Reorder Number	Device Diameter at Descending Aorta (mm) A	Device Diameter at Pulmonary Artery (mm) B	Retention Skirt Diameter (mm) C	Device Length (mm) D	Model/Reorder Number	Minimum Recommended Sheath Size	Sheath Inner Diameter (mm [inch])	Sheath Outer Diameter (mm [inch])
9-PDA-003	5	4	9	5	9-ITV05F180/60	5 F	1.83 [0.072]	2.51 [0.099]
9-PDA-004	6	4	10	7				
9-PDA-005	8	6	12	7	9-ITV06F180/80	6 F	2.11 [0.083]	2.79 [0.110]
9-PDA-006	10	8	16	8				
9-PDA-007	12	10	18	8				
9-PDA-008	14	12	20	8	9-ITV07F180/80	7 F	2.44 [0.096]	3.18 [0.125]
9-PDA-009	16	14	22	8				




## DEVICE SPECIFICATIONS AND DELIVERY SYSTEM COMPATIBILITY

### Amplatzer™ Duct Occluder II

### Amplatzer™ TorqVue™ LP Delivery System

Model/Reorder Number	Waist Diameter (mm) A	Device Length (mm) B	Disc Diameter (mm) C	Model/Reorder Number	Minimum Recommended Catheter Size	Catheter Inner Diameter (mm [inch])	Catheter Outer Diameter (mm [inch])
9-PDA2-03-04	3	4	9				
9-PDA2-03-06	3	6	9	9-TVLP4F90/060 or 9-TVLP4F90/080	4 F	1.17 [0.046]	1.40 [0.055]
9-PDA2-04-04	4	4	10				
9-PDA2-04-06	4	6	10				
9-PDA2-05-04	5	4	11				
9-PDA2-05-06	5	6	11	9-TVLP5F90/060 or 9-TVLP5F90/080	5 F	1.50 [0.059]	1.73 [0.068]
9-PDA2-06-04	6	4	12				
9-PDA2-06-06	6	6	12				



## DEVICE SPECIFICATIONS AND DELIVERY SYSTEM COMPATIBILITY

### Amplatzer™ Piccolo™ Occluder

### Amplatzer™ TorqVue™ LP Catheter

Model/Reorder Number	Disc Diameter (mm) A	Waist Diameter (mm) B	Length between discs (mm) C	Model/Reorder Number	Minimum Recommended Catheter Size	Catheter Inner Diameter (mm [inch])	Catheter Outer Diameter (mm [inch])
9-PDAP-03-02-L	4.00	3.00	2.00				
9-PDAP-03-04-L	4.00	3.00	4.00				
9-PDAP-03-06-L	4.00	3.00	6.00				
9-PDAP-04-02-L	5.25	4.00	2.00				
9-PDAP-04-04-L	5.25	4.00	4.00	9-TVLP4F90/080	4 F	1.17 [0.046]	1.40 [0.055]
9-PDAP-04-06-L	5.25	4.00	6.00				
9-PDAP-05-02-L	6.50	5.00	2.00				
9-PDAP-05-04-L	6.50	5.00	4.00				
9-PDAP-05-06-L	6.50	5.00	6.00				

## ANCILLARY PRODUCT SPECIFICATIONS

### Amplatzer™ Guidewire

Model/Reorder Number	Diameter (inch)	Body	Tip Description	Useable Length (cm)
9-GW-001	0.035	Super Stiff	7.5 mm, Modified J-tip	260

## LATEX-FREE INFORMATION

These Amplatzer™ products do not contain latex.

**CAUTION:** This product is intended for use by or under the direction of a physician. Prior to use, reference the Instructions for Use, inside the product carton (when available) or at [eifu.abbottvascular.com](http://eifu.abbottvascular.com) or at [medical.abbott/manuals](http://medical.abbott/manuals) for more detailed information on Indications, Contraindications, Warnings, Precautions and Adverse Events.

Information contained herein for **DISTRIBUTION in Australia and New Zealand ONLY**. Always check the regulatory status of the device in your region.

Illustrations are artist's representations only and should not be considered as engineering drawings or photographs. Photo(s) on file at Abbott.

### Abbott

Abbott Medical Australia Pty Ltd, 299 Lane Cove Road, Macquarie Park, NSW 2113, Ph: 1800 839 259  
Abbott Medical New Zealand Ltd, 4 Pacific Rise, Mount Wellington, Auckland 1060, Ph: 0800 756 269

™ Indicates a trademark of the Abbott group of companies.

‡ Indicates a third-party trademark, which is property of its respective owner.

© 2023 Abbott. All Rights Reserved. MAT-2304926 v1.0 | Item approved for ANZ use only.

