

# NAVITOR™ TAVI SYSTEM

SMART SEALING. EXCEPTIONAL STABILITY. UNCOMPROMISED ACCESS.



## Ordering Information

Navitor™ Transcatheter Aortic Valve

CATALOG NUMBER	VALVE SIZE (MM)	ANNULUS USE RANGE DIAMETER (MM) <sup>1</sup>	ANNULUS AREA (MM <sup>2</sup> ) <sup>1</sup>	ANNULUS PERIMETER (MM) <sup>1</sup>
NVTR-23	23	19–21	277–346	60–66
NVTR-25	25	21–23	338–415	66–73
NVTR-27	27	23–25	405–491	72–79
NVTR-29	29	25–27	479–573	79–85

FlexNav™ Delivery System

CATALOG NUMBER	EQUIVALENT INTEGRATED SHEATH DIAMETER (F) <sup>1</sup>	OUTER DIAMETER (MM) <sup>1</sup>	INTEGRATED SHEATH WORKING LENGTH (CM) <sup>1</sup>	WORKING LENGTH (CM) <sup>1</sup>	VASCULAR ACCESS DIAMETER (MM) <sup>1</sup>
FNAV-DS-SM	14	6.0	30	107	≥ 5.0
FNAV-DS-LG	15	6.3	30	107	≥ 5.5

Navitor™ Loading System

CATALOG NUMBER	
NVTR-LS-SM	The Navitor™ loading system facilitates valve preparation/loading onto the FlexNav™ delivery system. The loading system includes a loading funnel, loading base, base insert, loading tube, stent guide and leaflet tester. The LG loading system contains an extra base insert to facilitate loading of the 29 mm valve. All base inserts are labeled with the corresponding valve size. In addition, the base insert and loading tube are color-coded to differentiate between the different configurations.
NVTR-LS-LG	

1. Abbott data on file 90734545/90465559.

CAUTION: This product is intended for use by or under the direction of a physician. Prior to use, reference the Instructions for Use, inside the product carton (when available) or at [eifu.abbottvascular.com](http://eifu.abbottvascular.com) or at [medical.abbott/manuals](http://medical.abbott/manuals) for more detailed information on Indications, Contraindications, Warnings, Precautions and Adverse Events.

Information contained herein for **DISTRIBUTION in New Zealand ONLY**. Illustrations are artist's representations only and should not be considered as engineering drawings or photographs. Photo(s) on file at Abbott.

Abbott Medical New Zealand Ltd, Ground Floor, Bldg D, 4 Pacific Rise, Mount Wellington, Auckland 1060 Ph: 0800 656 233  
W: [www.aus.abbott](http://www.aus.abbott) E: [ANZStructuralHeart@abbott.com](mailto:ANZStructuralHeart@abbott.com)

™ Indicates a trademark of the Abbott group of companies.

‡ Indicates a third-party trademark, which is property of its respective owner.

©2021 Abbott. All rights reserved.

MAT-2108331 v1.0 | Item approved for NZ use only.

