




# Quartet<sup>TM</sup> Family of LV Leads



# PROVEN QUARTET™ LV LEAD PERFORMANCE

with the most Quadripolar lead options to match a patient's anatomy.

## BUILT ON THE PROVEN QUARTET™ LV LEAD PLATFORM AND FEATURING:

- Four uniquely spaced electrodes to provide flexibility in pacing location
- Small 4.7 F lead body and 4.0 F lead tip
- Optim™ lead insulation




The Quartet™ Family of LV leads offers more distal shape options including the Large-S, Small-S and the Double Bend and more total electrode spacing options including 40, 47 and 60 mm to offer not only successful CRT Implant but true quadripolar flexibility and options to help patient optimization.




**1457Q Double-Bend Shape  
Quartet™ LV Lead**

- May be ideal for medium size and length target vessels
- Total spacing 47 mm (20-10-17)




**1456Q Small-S Shape  
Quartet™ LV Lead**

- May be ideal for shorter target vessels
- Total spacing 40 mm (20-10-10)




**1458QL Traditional S Shape (wide spacing)  
Quartet™ LV Lead**

- May be ideal for larger size and longer length target vessels
- Total spacing 60 mm (20-27-13)



**1458Q Traditional S Shape  
Quartet™ LV Lead**

- Superb deliverability with exceptional stability and performance
- Total spacing 47 mm (20-10-17)



For more information regarding specific product combinations and associated MRI scan parameters, please see the MRI Ready Systems Manual at [manuals.sjm.com](http://manuals.sjm.com).

**Abbott**

One St. Jude Medical Dr., St. Paul, MN 55117 USA, Tel: 1 651 756 2000

[SJM.com](http://SJM.com)

St. Jude Medical is now Abbott.

**Rx Only**

**Brief Summary:** Prior to using these devices, please review the Instructions for Use for a complete listing of indications, contraindications, warnings, precautions, potential adverse events and directions for use.

<sup>™</sup> Indicates a trademark of the Abbott group of companies.

‡ Indicates a third party trademark, which is property of its respective owner.

© 2018 Abbott. All Rights Reserved.

25475-SJM-QRT-0116-0019a(1) | Item approved for U.S. use only.

